

Ανυψωτικό μασάζ προσώπου με RF ANLAN 01-ADRY51-02A

ANLAN · EAN 6977254876078

Machine translation from Polski. In case of doubt the manufacturer's original document prevails.

Lifting facial massager with RF ANLAN 01-ADRY51-02A

User Manual

This document is a translation of the original user manual created by the manufacturer.

Detailed information about the distributor/manufacturer warranty conditions is available on the website <https://serwis.innpro.pl/gwarancja>.

The product should be regularly maintained (cleaned) by the user or through specialized service points at the user's cost and responsibility. In the absence of information about necessary periodic or service maintenance actions in the user manual, the user should regularly, at least once a week, assess the physical condition of the product against that of a physically new product. If any deviation is detected or noted, urgent maintenance (cleaning) or service steps should be taken. Failure to properly maintain (clean) and respond upon detecting a state of deviation may lead to permanent damage to the product. The manufacturer is not responsible for damages resulting from neglect.

Warnings and Safety Information

All information regarding the use of the product is contained in the user manual. Before you start using it, familiarize yourself with its content and follow the instructions provided. Also, review the following information before use:

Usage Warnings

Risk of product damage:

- Use the devices only for their intended purpose, avoiding use on the body in a manner inconsistent with the instructions.
- Non-waterproof devices should not be submerged in water or used in high humidity conditions.
- Do not expose products to extreme temperatures or leave them in direct sunlight for extended periods.
- For battery-operated devices, ensure that charging is done only with the supplied cable or a compatible adapter.

User hazards:

- Avoid using devices on irritated, damaged skin, in the case of open wounds, active inflammatory conditions, or allergic reactions.
- Do not use devices that emit microcurrents, ultrasound, or heat if you have a pacemaker, metal implants, or are pregnant.
- If you experience an allergy, discomfort, or other concerning symptoms, immediately stop using the device and consult a dermatologist.

Information on Proper Use

Recommendations for use:

- Before using devices such as massagers or blackhead vacuums, thoroughly clean your face and apply appropriate cosmetics (e.g., gel or water) to assist the treatment.
- Deep cleansing devices, such as cavitation peels, should be used a maximum of 2-3 times a week to avoid excessive skin irritation.
- Use the devices according to the manufacturer's recommendations, not exceeding the recommended operating time on one area of the face (usually 1-2 minutes).
- Choose the appropriate intensity level, especially if you have sensitive skin.

Maintenance and Storage:

- After each use, clean the device with a damp cloth or another tool in accordance with the manufacturer's recommendations. For waterproof devices, you can rinse them under running water.
- Regularly check the condition of accessories, such as peeling tips or brushes, and replace them according to the manufacturer's recommendations.
- Store devices in a dry, clean place, avoiding contact with moisture and dust.

Additional Precautions

User Adaptation:

- Choose a device suitable for your skin type and individual needs. Features such as heating or microcurrent massage may be contraindicated for delicate or capillary skin.
- The use of devices by individuals with dermatological issues, such as rosacea, should be done after consulting a physician.

Damage and Service:

- If you notice any damage to the device (e.g., cracked casing, issues with turning on, unstable operation), immediately stop using it and contact an authorized service center.
- Use only original components supplied by the manufacturer for part and accessory replacements.
- Regularly check the operation of battery-operated devices and charge them according to

the recommendations to avoid performance drop or risk of damage.

Declaration of Conformity

The products have been designed in accordance with the safety requirements specified in Regulation (EU) 2023/988. They meet the relevant technical standards, including those concerning user safety, electromagnetic compatibility, and material durability.

Safety Measures

Before starting to charge, check that the device contacts are clean. Never leave the device unattended during use and charging. Ensure that you can quickly disconnect the device from the power source in an emergency.

Never expose the device to high temperatures. Charge the device in a dry and well-ventilated area away from flammable materials, maintaining a minimum distance of 1m from other objects. Never cover the device while charging. Never use a power supply, charging station, cables, etc., without the manufacturer's recommendation and certification. Take care of your property; the device is equipped with cells that are difficult to extinguish, so equip yourself with a fire blanket.

LI-ION Battery The device is equipped with a LI ION (lithium-ion) battery, which, due to its physical and chemical structure, ages over time and with use. The manufacturer specifies the maximum operating time of the device under laboratory conditions, where optimal working conditions for the device exist, and the battery itself is new and fully charged. The actual operating time may differ from the declared one in the offer and this is not a defect of the device but a feature of the product. To maintain the maximum lifespan of the battery, it is not recommended to discharge it below 3.18V or 15% of the total capacity. Lower values, such as 2.5V for the cell, permanently damage it and this is not covered by the warranty. If the battery or the entire device is not used for more than one month, the battery should be charged to 50% and its charge level should be checked periodically every two months. Store the battery and device in a dry place, away from sunlight and negative temperatures.

LI-PO Battery The device is equipped with a LI PO (lithium-polymer) battery, which, due to its physical and chemical structure, ages over time and with use. The manufacturer specifies the maximum operating time of the device under laboratory conditions, where optimal working conditions for the device exist, and the battery itself is new and fully charged. The actual operating time may differ from the declared one in the offer and this is not a defect of the device but a feature of the product. To maintain the maximum lifespan of the battery, it is not recommended to discharge it below 3.5V or 5% of the total capacity. Lower values, such as 3.2V for the cell, permanently damage it and this is not covered by the warranty. If the battery or the entire device is not used for more than one

month, the battery should be charged to 50% and its charge level should be checked periodically every two months. Store the battery and device in a dry place, away from sunlight and negative temperatures. Environmental Protection Used electronic equipment marked in accordance with the European Union directive must not be disposed of with other municipal waste. It is subject to selective collection and recycling at designated points. By ensuring its proper disposal, you prevent potential negative consequences for the natural environment and human health. The system for collecting used equipment complies with locally applicable environmental protection regulations regarding waste disposal. Detailed information on this topic can be obtained from the municipal office, waste management facility, or the store where the product was purchased. The product meets the requirements of the so-called New Approach directives of the European Union (EU) regarding issues related to safety of use, health protection, and environmental protection, specifying hazards that should be detected and eliminated. Importer: INNPRO Robert Błędowski sp. z o.o. Rudzka 65c 44-200 Rybnik, Poland tel. +48 533 234 303 hurt@innpro.pl www.innpro.pl Manufacturer: Shenzhen ANLAN Life Technology Co., Ltd 2204.Bldg.4, Cloud Park Phase 2, BanTian St. / 518129 Longgang, Shenzhen, Guangdong, China support@anlan.com Representative in the EU: LTI Testing & Certification France S.A.S. Rue de Rome 154 13006 Marseille, France info@lti-lab.com Product battery specifications Feature Value Battery category Portable battery Battery type Single Cell Li-ion Battery Contains Cadmium/Lead (>0.002% Cd, >0.004% Pb) No Net weight (kg) 0.0198 Capacity (mAh) Battery Capacity: 1100mAh