

Αποτριχωτική μηχανή IPL ANLAN 02-ATMY52-ORE

ANLAN · EAN 6977254876764

Machine translation from Polski. In case of doubt the manufacturer's original document prevails.

This document is a translation of the original user manual created by the manufacturer. Detailed information about the distributor/manufacturer warranty conditions is available on the website <https://serwis.innpro.pl/gwarancja>. The product should be regularly maintained (cleaned) at the user's expense and responsibility, either by the user themselves or through specialized service points. In the absence of information about necessary periodic or service maintenance actions in the user manual, the user should regularly, at least once a week, assess the physical condition of the product against that of a physically new product. If any deviation is detected or noted, urgent maintenance (cleaning) or service steps should be taken. Failure to properly maintain (clean) and respond upon detecting a state of deviation may lead to permanent damage to the product. The warranty provider is not responsible for damages resulting from neglect.

Warnings and safety information

All information regarding the use of the product can be found in the user manual. Before you start using it, familiarize yourself with its content and follow the guidelines contained within. Also, review the following information before use:

Warnings regarding usage

Risk of product damage:

- Products should only be used for their intended purpose.
- Do not use devices on damaged, irritated, or burned skin.
- In the case of IPL depilators, avoid using them on tanned skin to prevent discoloration and irritation.
- Electric brushes and graters should be regularly cleaned and checked for mechanical damage, such as cracks or worn heads.
- Do not immerse electrical devices in water unless the manufacturer has marked them as waterproof.
- Do not leave devices in a humid environment unless they are intended for such use.

User hazards:

- IPL depilators may cause skin irritation if used improperly. Conduct a patch test on a small

area of skin before first use.

- EMS devices (electric massagers) are not suitable for individuals with pacemakers, pregnant women, or people with skin diseases.
- Electric foot graters should only be used on dry skin. Avoid contact with open wounds.
- Avoid using heating devices on sensitive areas of the body to prevent burns.

Risks associated with electronic components:

- Use only original chargers and accessories provided by the manufacturer.
- Before each use, ensure that power cords and batteries are in good technical condition.
- If you notice damage to the cords or battery, immediately cease usage.

Information on proper usage

Recommendations for use:

- IPL depilators should be used on clean, dry skin, avoiding sensitive areas such as around the eyes. Regularity of treatments ensures the best results.
- Clean the massage brush heads after each use to prevent the buildup of bacteria. Adjust the intensity of the massage to the individual needs of the skin.
- Use interchangeable abrasive elements of graters and foot files according to the manufacturer's recommendations. After each use, clean the device of dead skin residue.
- For EMS devices, place electrodes only on clean, dry skin, avoiding sensitive areas and open wounds.

Maintenance and storage

- Regularly clean devices after each use using a damp cloth or dedicated cleaning agents.
- Store devices in a dry place, avoiding extreme temperatures and direct sunlight.
- Turn off battery-operated or rechargeable devices after use and charge them according to the user manual.

Additional precautions

User adaptation:

- Regularly monitor skin reactions to the use of devices. In case of irritation or redness, discontinue use and consult a dermatologist.
- Choose intensity settings and operating modes appropriate for skin type and sensitivity level.
- Before using heating devices or those operating at high intensity, consult a doctor, especially in the case of skin conditions or health issues.

Damage and service:

- If you notice leaks, operational issues, electrical discharges, or other damages, immediately cease usage and contact an authorized service center.
- For repairs and parts replacement, use only original components approved by the manufacturer.
- Do not make unauthorized modifications to devices, as this may affect their safety and void the warranty.

Declaration of conformity

The products have been designed in accordance with the safety requirements specified in Regulation (EU) 2023/988.

They meet the relevant technical standards, including those related to user safety, electromagnetic compatibility, and material durability. Environmental Protection Used electronic equipment marked in accordance with the European Union directive must not be disposed of together with other municipal waste. It is subject to selective collection and recycling at designated points. By ensuring its proper disposal, you prevent potential negative consequences for the natural environment and human health. The system for collecting used equipment complies with local environmental protection regulations regarding waste disposal. Detailed information on this topic can be obtained from the municipal office, cleaning company, or the store where the product was purchased. The product meets the requirements of the so-called New Approach directives of the European Union (EU) concerning issues related to user safety, health protection, and environmental protection, which define hazards that should be detected and eliminated. Importer: INNPRO Robert Błędowski sp. z o.o. Rudzka 65c 44-200 Rybnik, Poland tel. +48 533 234 303 hurt@innpro.pl www.innpro.pl Manufacturer: Shenzhen ANLAN Life Technology Co., Ltd 2204.Bldg.4, Cloud Park Phase 2, BanTian St. / 518129 Longgang, Shenzhen, Guangdong, China support@anlan.com EU Representative: LTI Testing & Certification France S.A.S. Rue de Rome 154 13006 Marseille, France info@lti-lab.com