


ANLAN

MIRACLE 4-IN-1 MASSAGER GUA SHA BEAUTY SET

USER MANUAL



Thank you for purchasing the ANLAN Miracle Beauty Gua Sha Set. Before using this product, please read this user manual carefully for future reference. In addition, this manual also serves as the warranty, so please keep it properly.

Safety Precautions

[Contraindications]

- Please read the instructions carefully. Do not use this product on the same area continuously for extended periods of time to avoid skin damage, frostbite, burns, or muscle pain caused by excessive massage of the same area.
- Do not use this product while performing delicate operations or while driving.
- Do not use this product if the massage area has undergone surgery or if there are any skin injuries, redness, swelling, inflammation, bruising, or other discomfort.
- Do not use this product if you are sensitive to metal or heat, or if you have a history of heart disease, acute illness, skin disease, malignant tumor, bleeding disorder, osteoporosis, soft tissue disease, or stroke.
- Do not use this product if you have motor impairment, physical or mental impairment, or intellectual disability. If you need to use this product, please consult a doctor.
- It is not recommended to use this product during pregnancy or menstruation. If you need to use this product, please consult a doctor.
- Store this product out of the reach of children. 8. Not for use by children under six years old; children aged 6-12 must use under the guidance and supervision of a parent or other competent adult.

[Warm Tips]

- This product is a home electronic health product and is not intended to replace a medical device. It is not intended for diagnostic, therapeutic, or other

commercial purposes.

- Before use, please remove any jewelry or other items that may interfere with the massage area.
- This product's massage head has a metal surface. Users with allergies to such surfaces are prohibited from using it.
- This product has both heating and cooling surfaces. It should not be used by individuals insensitive to heat and cold, infants, or individuals unable to care for themselves.
- When using this product, pay attention to the surface temperature of the heated metal and ensure that your skin can tolerate it before use to avoid burns.
- If any odor remains after using this product with skincare products such as essential oils or creams, wipe it clean with a clean, damp cloth to keep it clean. Disinfect with 75% alcohol before using it on others to prevent infection.
- Do not use chemical solvents to clean the body, power supply, or button cover of this product.
- Discontinue use if you experience any pain or discomfort in the massage area.
- If any abnormality occurs during use, immediately turn off the product, disconnect the power cord, and contact customer service. Do not attempt to disassemble or repair the product yourself.
- When using the product in other specific environments, please comply with safety regulations.
- Please charge the product with a power adapter that matches the product's input voltage and has a maximum output power of 15W. Pay attention to electrical safety. It is recommended not to use the product with other charging devices.

- Do not disassemble the product yourself. Any damage caused by any operation performed by anyone other than our technicians will be the sole responsibility of the user.
- Keep this product away from water, fire, and corrosive environments.
- Excessive use may cause physical discomfort or skin damage. Therefore, upon first use, observe your body's reactions and skin changes. We recommend shortening the use time.

[Before Use]

- When using this product, maintain a comfortable position.
- Before use, remove any jewelry, etc. to avoid discomfort during massage.
- If you experience any skin pain or discomfort, discontinue use and seek medical attention immediately.

[After Use]

- The massaged area may briefly become red and warm. This is normal and will resolve after a while.


[Cleaning and Protection]

- To clean the exterior of the product, wipe it with a cloth dampened with water and wrung out thoroughly. Do not rinse the body directly with water. Rinse the magnetic replacement head with clean water and then dry it (do not soak it in water).
- ※ Keep this product away from corrosive environments.

[Other]

- During normal use, the product may make a buzzing sound. This is part of normal operation.

Introduction to product functional



Specifications


Name	Beauty Device Comb Massager 4-in-1 Set
Model	DYM-072
Dimensions	Main unit (without head): Appr 131*39*30mm Metal Mountain Massage Head: Appr 101*41*27mm Wide-tooth Comb Head: Appr 101*41*25mm Close-tooth Comb Head: Appr 101*41*31mm 7D Metal Massage Ball: Appr 101*41*22mm

Net Weight	appr 387g
Battery capacity	1200mAh
Input	5V=1A
Battery life	Metal Mountain Massage Head Cool: Appr 90 minutes; Other: 180 minutes
Charge time	about 120 minutes
Material	ABS+PC+zinc alloy+ stainless steel
Rated input power	5W

Packing List



Main unit*1



magnetic replaceable massage heads*4



type-c charging cable*1



instruction manual*1




velvet bag *1

How to use the product

Before use, please remove any jewelry or accessories from the massage area to avoid discomfort. This

product cannot be used while charging. Please unplug the power adapter after charging is complete.

1.Metal Mountain Massage Head (Hot/Cold)



① Power On

Press and hold the button for 0.5 seconds to power on the device. If no massage head is attached, one blue light will illuminate. After one minute of standby, the device will shut down after two beeps.

② Mode Switching

[Cold Compress]

This head defaults to cold compress mode upon powering on. Setting the first position, one blue light will illuminate, and the cold compress power will be 70%. Short-pressing the power button will shift the device to the second position, with two blue lights illuminating and the cold compress power at 85%. Short-pressing the power button will shift the device to the third position, with three blue lights illuminating and the cold compress power at 100%. (At room temperature of 26°C, setting the first position is 22°C ± 2°C; the second position is 20°C ± 2°C; and the third position is 18°C ± 2°C. Cooling performance may vary depending on ambient temperature.)

[Hot Compress]

When in the third position of the cold compress mode, short-pressing the power button will shift the device to hot compress mode. One red light will illuminate, and the first position of the hot compress mode will be activated. The massage head temperature rises to approximately $38^{\circ}\text{C} \pm 2^{\circ}\text{C}$. Short-press the power button to shift to level 2. Two red lights illuminate, and the massage head temperature rises to approximately $42^{\circ}\text{C} \pm 2^{\circ}\text{C}$. Short-press the power button to shift to level 3. Three red lights illuminate, and the massage head temperature rises to approximately $46^{\circ}\text{C} \pm 2^{\circ}\text{C}$. Short-press the power button again to switch to cold compress mode. Short-press the button to cycle through cold compress levels 1-3 and hot compress levels 1-3.

③ If the massage head is removed, the buzzer will beep three times and the massage head will enter standby mode.

※ CAUTION: Avoid using hot or cold compresses on the same area for extended periods to avoid burns and frostbite.

④ Timer Function

Hot compress mode automatically times out for 10 minutes, followed by a beep. Cold compress mode only times out for 3 minutes, followed by a beep.

⑤ Shutdown Press and hold the button for 0.5 seconds to power off the product.

⑥ Usage Tips


[Recommended usage time] Approximately 3–8 minutes

[Recommended frequency of use] Approximately 3–5 times per week

[Suitable areas] Face, neck, shoulders, back, and limbs Can be used with massage oil or body lotion to enhance the smoothness and comfort of scraping.

(Effectiveness varies from person to person; the above information is for reference only.)

2. Wide-tooth Comb Head (purple light + microcurrent)



① Power On

Press and hold the button for about 0.5 seconds to power on the product. If the massage head is not attached, one blue light will illuminate. After one minute of standby, the product will shut down after two beeps.

② Mode Switching

[Microcurrent Mode]

With the massage head attached, power on with one purple light illuminated, indicating the product is in the first position (65% brightness). Short press the button to activate the product, with two purple lights illuminating, indicating the product is in the second position (85% brightness). Short press the button to activate the product,

with three purple lights illuminating, indicating the product is in the third position (100% brightness). Short press the power button to cycle through the power modes.

③ If the massage head is removed, the buzzer will beep three times shortly, and the product will enter standby mode.

④ **Timer**

This mode automatically sets a 10-minute timer, followed by a beep to shut down the product. Short press the button for about 0.5 seconds to power off the product.

⑤ **Power Off**

Press and hold the button for about 0.5 seconds to power off the product.

⑥ **Usage Tips**

【Recommended usage time】 Approximately 5–8 minutes


【Recommended frequency of use】 Approximately 3–5 times per week

【Applicable areas】 Scalp, neck, shoulders, and back

【Directions】 Recommended for use with conductive gel or scalp treatment serum to enhance the microcurrent delivery effect.

(Effectiveness varies from person to person; the following is for reference only.)

3. Close-tooth Comb Head (Red Light, Purple Light + Microcurrent)



① **Power On**

Press and hold the button for about 0.5 seconds to power on the device. If the massage head is not attached, one blue light will illuminate. The device will then standby for one minute and shut down after two beeps.

② **Mode Switching**

【Red Light Mode】

With this head attached, power on, one red light will illuminate, indicating the device is in the first position (65% brightness). Short press the button, two red lights will illuminate, indicating the device is in the second position (85% brightness). Short press the button, three red lights will illuminate, indicating the device is in the third position (100% brightness). Short press the power button to switch to Purple Light + Microcurrent Mode.

【Purple Light + Microcurrent Mode】

With this head attached, power on, one purple light will illuminate, indicating the device is in the first position and 65% brightness. Short press the button, two purple lights will illuminate, indicating

the device is in the second position and 85% brightness. Short press the button, three purple lights will illuminate, indicating the device is in the third position and 100% brightness. Short press the power button to cycle through the modes and modes.

③ If you remove the massage head, the buzzer will beep three times shortly and the massager will enter standby mode.

④ Timer Function

This mode sets the timer to 10 minutes, followed by a single beep to shut down the massager.

⑤ Power Off

Press and hold the button for about 0.5 seconds to shut down the massager.

⑥ Usage Tips

【Recommended usage time】 Approximately 5–8 minutes per session.


【Recommended frequency of use】 Approximately 3–5 times per week.

【Applicable areas】 Scalp and hairline.

【How to use】 Can be used with hair growth serums, scalp oils, etc. to promote nutrient absorption.

(Effectiveness varies from person to person; the following is for reference only.)

4.7D Metal Massage Ball (Purple Light + Microcurrent)



① Power On

Press and hold the button for about 0.5 seconds to power on the product. If the massage head is not attached, one blue light will illuminate. After one minute of standby, the product will shut down after two beeps.

② Mode Switching

【Microcurrent Mode】

With the massage head attached, one purple light will illuminate, indicating the first setting and 65% brightness. Short-press the button to illuminate two purple lights, indicating the second setting and 85% brightness. Short-press the button to illuminate three purple lights, indicating the third setting and 100% brightness. Short-press the power button to cycle through the settings.

③ If you remove the massage head, the buzzer will beep three times shortly and the massager will enter standby mode.

④ Timer Function

This mode automatically sets the timer to 10 minutes, after which it will shut down with a single beep.

⑤ Power Off

Press and hold the button for about 0.5 seconds to turn the massager off.

⑥ Usage Tips

[Recommended usage time] Approximately 5–8 minutes

[Recommended frequency of use] Approximately 3–5 times per week

[Applicable areas] Facial contours, jawline, décolleté

[How to use] It is recommended to use with conductive gel or facial serum for enhanced conductivity and comfort.

(Effectiveness may vary from person to person; the following is for reference only.)

Recycling and Disposal Methods



Do not disassemble unless authorized by a qualified professional.

【 I. Battery Removal Procedure 】

1. Disconnect the external power supply and remove the product casing.
2. Unplug the power cord from the battery.
3. Remove the battery.
4. Dispose of used batteries at a designated recycling point.
5. Please comply with local environmental laws and disposal guidelines to protect the environment.

【 II. Storage 】

If the storage period exceeds three months, unpack and inspect the device, recharge it, and then use it again. Avoid placing it against the floor, walls, or ceiling. Ensure the room is well ventilated and protected from strong sunlight and corrosive gases.

【 III. Transportation 】

The product packaging is equipped with simple shockproofing for transportation by air, rail, road, and ship. Protect it from rain, snow, inversion, and impact.

【 IV. Power Supply Inspection 】

Check the power supply. Do not operate the device if the power voltage is outside the specified range.

【 V. Maintenance 】

Take measures to prevent moisture, rodents, and contamination. Clean the device with the power off. After use, store the device in a safe place. Do not place heavy objects on top of it. In the event of a malfunction, consult a qualified technician for repair. Non-professionals are not allowed to disassemble or assemble the product.

【 VI. Disposal and Recycling 】

This product should not be disposed of with household waste. When you decide to dispose of this product, please comply with local environmental laws and guidelines.



Note

Before disposing of this product, the battery must be removed as required by law.

When removing the battery, the product must be powered off.

The battery must be disposed of safely and in accordance with laws and regulations.