


ANLAN

2-in-1 Guasha Beauty Device

USER MANUAL



※ Thank you for purchasing the ANLAN 2-in-1 Guasha Beauty Device. Before using this product, please read this user manual carefully for future reference. In addition, this manual also serves as a warranty, please keep it properly.

Safety Precautions

PROHIBIT

- Please read the instructions carefully. Do not use this product on the same area continuously for extended periods of time to avoid skin damage, frostbite, burns, or muscle pain caused by excessive massage of the same area.
- Do not use this product while performing delicate operations or while driving.
- Do not use this product if the massage area has undergone surgery or if you have skin injuries, redness, swelling, inflammation, bruising, or other discomfort.
- Do not use this product if you are sensitive to metal, Bianstone, or heat, or if you have a history of heart disease, acute illness, skin disease, malignant tumor, bleeding disorder, osteoporosis, soft tissue disease, or stroke.
- Do not use this product if you have motor delays, physical or mental impairments, or intellectual disabilities. If you need to use this product, please consult a doctor.
- It is not recommended to use this massager during pregnancy or menstruation. If you need to use this product, please consult a doctor.
- Store this product out of the reach of children. 8. Children under six years old are prohibited from using this product. Children aged 6-12 must use this product under the guidance and supervision of their parents or other adults with supervisory capabilities.


PROHIBIT

- This product is a home electronic health product and is not intended for use as a substitute for medical equipment. It is not intended for diagnosis, treatment, or other commercial purposes.
- Before use, please remove any jewelry or other items that may interfere with the massage area.
- This product's massage head has metal or Bianstone surfaces. Users with allergies to these materials are prohibited from use.
- This product has a heating surface and should not be used by individuals who are insensitive to heat, infants, or individuals who are unable to care for themselves.
- When using this product, pay attention to the surface temperature of the heated metal and Bianstone surfaces and ensure your skin can tolerate it before use to avoid burns.
- If any odor remains after using this product with skincare products such as essential oils or creams, wipe it clean with a clean, damp cloth to keep it clean. Disinfect with 75% alcohol before using it on others to prevent infection.
- Do not use chemical solvents to clean the body, power supply, or button cover of this product.
- Discontinue use if you experience any pain or discomfort in the massage area.
- If the product exhibits any abnormal behavior during use, immediately turn it off, disconnect the power cord, and contact customer service. Do not attempt to disassemble or repair the product yourself.
- When using the product in other specific environments, please comply with safety regulations.


⊘ PROHIBIT

- Please charge the product with a power adapter that matches the product's input voltage and has a charging output power of no more than 15W. Pay attention to electrical safety. It is recommended not to use the product with other charging devices.
- Do not disassemble the product yourself. Any damage caused by unauthorized operation by non-technical personnel will be the sole responsibility of the user.
- Keep the product away from water, fire, and corrosive environments.
- Excessive use may cause physical discomfort or skin damage. Therefore, during initial use, observe your body's reactions and skin changes. Shorter usage times are recommended.

Product Introduction



- ① Mode/Vibration button
- ② Power/Level button
- ③ Indicator




- ④ Guasha massage head
- ⑤/⑥ Beauty light
- ⑦ Meridian comb massage head

Specifications

Name	2 in 1 Guasha Beauty Device
Model	DYM-076
Input	5V ---1A
Battery Capacity	800mAh
Working Time	60min~150min
Charging Time	About 2h
Material	ABS+PC+Bianstone+Stainless Steel
Power	5W
Net Weight	170g
Size	About 109.6*80.6*32.2mm

Packing List



Device × 1



Type-C
charging cable × 1



User
manual × 1

How To Use

About Modes

Head	Guasha Head	Meridian Comb
Recommended cosmetics	Aroma oil, body lotion	Gel, scalp care essence
Recommended frequency of use	3-5 times a week 5-8 minutes each time	
Recommended areas for use	Face, neck, shoulders, back, limbs	Scalp, face, neck, shoulders, back
Level & Mode	Level 1: About 38°C, red light on, weak vibration intensity	Level 1: Red light on, weak vibration intensity & EMS intensity
	Level 2: About 42°C, blue light on, medium vibration intensity	Level 2: Blue light on, medium vibration intensity & EMS intensity
	Level 3: About 46°C, purple light on, strong vibration intensity	Level 3: Purple light on, strong vibration intensity & EMS intensity
	All levels are heating functions, only the temperature is different.	All levels are EMS functions, only the EMS intensity differs.

- ※When using the guasha head, avoid applying heat to the same area for a long time to avoid burns.
- ※The area may become red or warm for a short time after guasha treatment, but this is normal and will heal over time.
- ※The specific results vary from person to person, so the above information is for reference only.
- ※You may hear a buzzing sound during normal operation; this is a normal sound.


Operating

【Power ON/OFF】

Press and hold the power button  for about 0.5s, the power-on buzzer will beep for 0.5s, and the device will enter the Guasha heat compress mode. The power-off buzzer will beep for 1s.

【Mode Switching】

After turning on the device, short press the mode button to cycle between the two modes (Guasha heat compress / meridian comb microcurrent)



※In each mode, you can turn vibration on/off by pressing and holding the Mode/Vibration button.

【Level Adjustment】

Short press the power button to cycle through the three levels.



【Timing】

The device works for 10 minutes in each mode/each level, and it will automatically shut down after a "beep" sound.

Recommended Usage

Scalp

Select your preferred mode and slide in the direction of the arrow.

※This usage only applies to the meridian



Face

Select your preferred mode and slide in one direction of the arrow.

※Avoid the eyeballs.



Neck

Select your preferred mode and slide in one direction of the arrow.

※Avoid the throat.




Arms

Select your preferred mode and slide in the direction of the arrow.



Calf

Select your preferred mode and slide in the direction of the arrow.



About Charging

Low Battery	The indicator light will flash red five times, a warning beep will sound five times, and then the power will automatically turn off.
Charging	The indicator light flashes red.
Fully charged	The indicator light will remain blue.

※The device cannot be powered on while charging. Additionally, if you plug in a USB cable to charge the device while it's running, it will switch to charging mode.

*Please use a power adapter that matches the product's input voltage and has a charging output power rating of no more than 15W to charge the device. Please pay attention to electrical safety. We recommend that you do not use this device with other charging devices.

Q & A

If a fault occurs during use, please check it according to the following methods. If the fault cannot be corrected or the product still does not operate as expected, please contact customer service for repair. Do not disassemble the product yourself!

Faults	Reasons	Solutions
Unable to turn on	Is the battery low?	If so, please charge it in time
The massage head is non-conductive	Strong personal resistance/dry climate	If so, moisturize your skin
The massage head does not heat	Is the hot compress button turned off?	If so, please turn on the hot compress

Product Disposal



Non-professionals please do not disassemble

【Battery removal steps】

1. Disconnect the external power supply and remove the product casing.
2. Remove the battery.
3. Dispose of used batteries at a designated recycling point.
4. Please comply with local environmental laws and disposal guidelines to protect the environment.

【Storage】

If the storage period exceeds 3 months, the device should be inspected and recharged before use. Do not place the device close to the ground, walls, or ceiling. The room should be well ventilated and protected from strong sunlight and corrosive gases.

【Transportation】

The product packaging is equipped with simple shockproof facilities and is suitable for transportation by air, rail, road and ship. It should be protected from rain, snow splashing, inversion and collision.

【Power supply check】

Check the power supply. If the power supply voltage is out of the specified range, do not turn on the machine.

【Maintainance】

Take measures to prevent moisture, rodents, and contamination. Always turn off the instrument when cleaning. After use, place the instrument in a safe place and do not place heavy objects on top of it. If a malfunction occurs, consult a qualified technician. Non-professionals are not allowed to disassemble or assemble the product.

【Disposal and recycling】

This product should not be disposed of with household waste. If you decide to dispose of this product, please follow local environmental laws and guidelines.

Note:

※Before discarding this product, the battery must be removed as required by law.

When removing the battery, the product must be powered off.

The battery must be disposed of safely and in accordance with local laws and regulations.